



# Rabbit Anti-TYLCV Polyclonal antibody (CABT-CS196)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	There is no cross reaction with healthy plant tissues such as corn and wheat.
Target	TYLCV
Immunogen	TYLCV
Source/Host	Rabbit
Species Reactivity	TYLCV
Conjugate	unconjugated
Applications	ELISA (cap)
Format	Liquid
Size	500 tests
Buffer	PBS
Storage	Store at 4°C.

## BACKGROUND

Keywords	TYLCV ; Tomato Yellow Leaf Curl Virus
----------	---------------------------------------