



Rabbit Anti-TYLCV Polyclonal antibody [AP] (CABT-CS197)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	There is no cross reaction with healthy plant tissues such as corn and wheat.
Target	TYLCV
Immunogen	TYLCV
Source/Host	Rabbit
Species Reactivity	TYLCV
Conjugate	AP
Applications	ELISA (det)
Format	Liquid
Size	500 tests
Buffer	PBS
Storage	Store at 4°C.

BACKGROUND

Keywords	TYLCV ; Tomato Yellow Leaf Curl Virus
----------	---------------------------------------