



Rabbit Anti-SPFMV Polyclonal antibody (CABT-CS202)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	There is no cross reaction with healthy plant tissues such as rice and other host plants.
Target	SPFMV
Immunogen	SPFMV
Source/Host	Rabbit
Species Reactivity	SPFMV
Conjugate	unconjugated
Format	Liquid
Size	500 tests
Buffer	PBS
Storage	Store at 4°C.

BACKGROUND

Keywords	SPFMV; Sweat Potato Feathery Mottle Virus
----------	---