



# Rabbit Anti-SPFMV Polyclonal antibody [AP] (CABT-CS203)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	There is no cross reaction with healthy plant tissues such as rice and other host plants.
Target	SPFMV
Immunogen	SPFMV
Source/Host	Rabbit
Species Reactivity	SPFMV
Conjugate	AP
Format	Liquid
Size	500 tests
Buffer	PBS
Storage	Store at 4°C.

## BACKGROUND

Keywords	SPFMV; Sweat Potato Feathery Mottle Virus
----------	---