



Mouse Anti-HPV 59 L1 monoclonal antibody (set) (CABT-CS280)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	The antibody set can be used for ELISA and screening appropriate mAbs to establish VLP vaccine in vitro specific activity assay.
Target	HPV 59 L1
Immunogen	E. coli expressed HPV59 L1 VLP
Isotype	IgG
Source/Host	Mouse
Species Reactivity	HPV 59
Clone	2E4, 4B6, 6A12, 6D9, 7D7, 8E1
Purification	Protein G purified
Conjugate	unconjugated
Applications	ELISA, Neut, Screening
Format	Liquid
Size	20 µg×6
Buffer	Tris-Gly
Preservative	None
Storage	Store at -20°C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Introduction	Human papillomavirus infection (HPV infection) is an infection caused by the human papillomavirus (HPV). About 90% of HPV infections cause no symptoms and resolve spontaneously within two years. However, in some cases, an HPV infection persists and results in either warts or precancerous lesions.
Keywords	HPV 59 L1; HPV; HPV 59; HPV L1