



Rabbit Anti-HPV 68-L1 Polyclonal antibody [HRP] (CABT-CS287)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody has different degrees of cross-reactivity to L1 proteins of other types of HPV
Target	HPV 68 L1
Immunogen	HPV 68 L1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HPV 68
Purification	Protein G purified
Conjugate	HRP
Applications	WB, ELISA
Format	Liquid
Concentration	1 mg/mL
Size	100 µg
Buffer	10mM PBS, 25% glycerol
Preservative	None
Storage	Store at -20°C. Avoid repeated freezing and thawing cycles.
Ship	Shipped with wet ice

BACKGROUND

Introduction	Human papillomavirus infection (HPV infection) is an infection caused by the human papillomavirus (HPV). About 90% of HPV infections cause no symptoms and resolve spontaneously within two years. However, in some cases, an HPV infection persists and results in either warts or precancerous lesions.
Keywords	HPV 68 L1; HPV; HPV 68; HPV L1