



# Mouse Anti-Gardnerella vaginalis Monoclonal antibody, clone 5082 (CABT-CS288)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	Gardnerella vaginalis
Target	Gardnerella vaginalis
Immunogen	Gardnerella vaginalis
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Gardnerella vaginalis
Purification	protein A chromatography
Conjugate	unconjugated
Applications	IF, ELISA, LFIA
Format	Liquid
Size	100 µg
Buffer	0.01 M phosphate buffered saline, pH 7.2
Preservative	0.1% sodium azide
Storage	This product contains no stabilizing proteins and should be stored at 2-8°C until ready for use.

## BACKGROUND

**Introduction**

Gardnerella vaginalis is a major species present in bacterial vaginosis. It is a small gram variable, anaerobic coccobacillus. Overgrowth of this bacterium occurs due to an imbalance of normal vaginal flora which can result from prescription antibiotics.

---

**Keywords**

Gardnerella vaginalis; G. vaginalis

---