



Goat Anti-S. cerevisiae Polyclonal antibody [Biotin] (CABT-CS295)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	S. cerevisiae
Target	S. cerevisiae
Isotype	IgG
Source/Host	Goat
Species Reactivity	S. cerevisiae
Purification	Affinity Purified
Conjugate	Biotin
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Store at 2-8°C

BACKGROUND

Keywords	Saccharomyces cerevisiae; S. cerevisiae
----------	---