



# Goat Anti-Human C-3 (B1C/B1A) Antiserum (CABT-CS346)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                     |
|--------------------|---------------------|
| Specificity        | Human C-3 (B1C/B1A) |
| Target             | C-3 (B1C/B1A)       |
| Source/Host        | Goat                |
| Species Reactivity | Human               |
| Conjugate          | unconjugated        |
| Applications       | ID                  |
| Format             | Liquid              |
| Size               | 0.5 lbs             |
| Buffer             | Antiserum           |

## BACKGROUND

|              |                                                                                                                                                                         |
|--------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Introduction | Complement component 3, often simply called C3, is a protein of the immune system. It plays a central role in the complement system and contributes to innate immunity. |
| Keywords     | Complement component 3; C3                                                                                                                                              |