



# Goat Anti-Human C-3 Inactivator (KAF) Antiserum (CABT-CS347)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	Human C-3 Inactivator (KAF)
Target	C-3 Inactivator (KAF)
Source/Host	Goat
Species Reactivity	Human
Conjugate	unconjugated
Applications	ID
Format	Liquid
Size	0.5 lbs
Buffer	Antiserum

## BACKGROUND

Introduction	Complement component 3, often simply called C3, is a protein of the immune system. It plays a central role in the complement system and contributes to innate immunity.
Keywords	Complement component 3; C3; C3 Inactivator; C3INA; Complement component 3 Inactivator