



Goat Anti-Human C-4 (B1E) Antiserum (CABT-CS349)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human C-4 (B1E)
Target	C-4 (B1E)
Source/Host	Goat
Species Reactivity	Human
Conjugate	unconjugated
Applications	ID
Format	Liquid
Size	0.5 lbs
Buffer	Antiserum

BACKGROUND

Introduction	Complement component 4 (C4), in humans, is a protein involved in the intricate complement system, originating from the human leukocyte antigen (HLA) system. It serves a number of critical functions in immunity, tolerance, and autoimmunity with the other numerous components. Furthermore, it is a crucial factor in connecting the recognition pathways of the overall system instigated by antibody-antigen (Ab-Ag) complexes to the other effector proteins of the innate immune response.
Keywords	Complement component 4; C4