



Goat Anti-Human C-5 (B1F) Antiserum (CABT-CS350)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human C-5 (B1F)
Target	C-5 (B1F)
Source/Host	Goat
Species Reactivity	Human
Conjugate	unconjugated
Applications	ID
Format	Liquid
Size	0.5 lbs
Buffer	Antiserum

BACKGROUND

Introduction	Complement component 5 is the fifth component of complement, which plays an important role in inflammatory and cell killing processes. This protein is composed of alpha and beta polypeptide chains that are linked by a disulfide bridge. An activation peptide, C5a, which is an anaphylatoxin that possesses potent spasmogenic and chemotactic activity, is derived from the alpha polypeptide via cleavage with a C5-convertase.
Keywords	Complement component 5; C5