



Goat Anti-Human C1q Antiserum (CABT-CS352)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human C1q
Target	C1q
Source/Host	Goat
Species Reactivity	Human
Conjugate	unconjugated
Applications	ID
Format	Liquid
Size	0.5 lbs
Buffer	Antiserum

BACKGROUND

Introduction	C1 is the first complement component in the cascade referred to as the classical pathway of complement. C1 binds to and is activated by antibodies bound to antigens (immune complexes) yielding a protease that initiates the cascade. C1 is actually a non-covalent complex of three different proteins (C1q, C1r and C1s) bound together in a calcium-dependent complex.
Keywords	C1 Complex; C1q; complement component 1q; C1