



Goat Anti-Human Properdin Antiserum (CABT-CS353)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human Properdin
Target	Properdin
Source/Host	Goat
Species Reactivity	Human
Conjugate	unconjugated
Applications	ID
Format	Liquid
Size	0.5 lbs
Buffer	Antiserum

BACKGROUND

Introduction	Properdin is plasma glycoprotein that activates the complement system of the innate immune system. This protein binds to bacterial cell walls and dying humans cells to stabilize the C3 and C5-convertase enzyme complexes to form an attack complex that lead to the lysis of the cell.
Keywords	complement factor properdin; CFP; complement factor P; properdin; BFD; PFC; PFD