



Human Anti-Human Mesothelin (amatuximab biosimilar) Monoclonal antibody, clone amatuximab (CABT-CS370)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	MSLN
Target	MSLN
Isotype	IgG1
Source/Host	Chimeric
Species Reactivity	Human
Clone	amatuximab
Purification	Purified from cell culture supernatant by affinity chromatography
Conjugate	unconjugated
Applications	ELISA
Reconstitution	Reconstitute with deionized water
Format	Powder
Size	50 µg, 100 µg
Buffer	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.
Preservative	0.1% Procline 300

Storage

Store at -20°C (Avoid repeated freezing and thawing)

BACKGROUND

Introduction

Mesothelin is a GPI-linked cell-surface glycoprotein, that is present on the surface of normal mesothelium and is overexpressed in many patients with epithelial ovarian cancer and malignant mesotheliomas.

Keywords

MSLN; Mesothelin; MPF

GENE INFORMATION

Entrez Gene ID

10232
