



# Human Anti-Human GPRC5D (Talquetamab biosimilar) Monoclonal antibody, clone Talquetamab (CABT-CS397)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	GPRC5D
<b>Target</b>	GPRC5D
<b>Isotype</b>	IgG4
<b>Source/Host</b>	Humanized
<b>Species Reactivity</b>	Human
<b>Clone</b>	Talquetamab
<b>Purification</b>	Purified from cell culture supernatant by affinity chromatography
<b>Conjugate</b>	unconjugated
<b>Applications</b>	FC
<b>Reconstitution</b>	Reconstitute with deionized water
<b>Format</b>	Powder
<b>Size</b>	50 µg, 100 µg
<b>Buffer</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.
<b>Preservative</b>	0.1% Procline 300

**Storage**

Store at -20°C (Avoid repeated freezing and thawing)

---

## BACKGROUND

**Introduction**

The protein encoded by this gene is a member of the G protein-coupled receptor family; however, the specific function of this gene has not yet been determined.

---

**Keywords**

GPRC5D

---

## GENE INFORMATION

**Entrez Gene ID**

55507

---