



Duck Anti-Pit viper venom Polyclonal Antibody (CABT-CS488)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Bothrops jararaca venom
Target	Bothrops jararaca venom
Isotype	IgY (Δ Fc)
Source/Host	Duck
Species Reactivity	B. jararaca
Purification	>95%
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Store at 2-8°C for short term. Store at -20°C for long term. Avoid repeated freezing-thawing cycles.

BACKGROUND

Introduction	The Brazilian pit viper or jararaca, Bothrops jararaca, recently renamed as Bothropoides jararaca by Fenwick et al., produces a highly toxic venom which effects hemostasis, causing consumptive coagulopathy and local and systemic hemorrhage
Keywords	Bothrops jararaca venom; Pit viper venom