



Duck Anti-EV71 Polyclonal Antibody (CABT-CS494)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	EV71
Target	EV71
Isotype	IgY (Δ Fc)
Source/Host	Duck
Species Reactivity	EV71
Purification	>99%
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Store at 2-8°C for short term. Store at -20°C for long term. Avoid repeated freezing-thawing cycles.

BACKGROUND

Introduction	Enterovirus A71 (EV-A71) is a significant human pathogen, especially in children. EV-A71 infection is one of the leading causes of hand, foot, and mouth diseases (HFMD), and can lead to neurological complications such as acute flaccid myelitis (AFM) in severe cases.
Keywords	EV71; EV 71; EV-71; Enterovirus 71