



Duck Anti-HBsAg Polyclonal Antibody (CABT-CS495)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	HBsAg
Target	HBsAg
Isotype	IgY (Δ Fc)
Source/Host	Duck
Species Reactivity	HBV
Purification	>95%
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Store at 2-8°C for short term. Store at -20°C for long term. Avoid repeated freezing-thawing cycles.

BACKGROUND

Introduction	HBsAg (also known as the Australia antigen) is the surface antigen of the hepatitis B virus (HBV). Its presence in blood indicates current hepatitis B infection.
Keywords	HBV Surface Antigen; hepatitis B Surface Antigen; HBsAg; HBV; Hepatitis B Virus; Major surface antigen