



Duck Anti-H. pylori Polyclonal Antibody (CABT-CS497)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	H. pylori
Target	H. pylori
Isotype	IgY (Δ Fc)
Source/Host	Duck
Species Reactivity	H. pylori
Purification	>99%
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Store at 2-8°C for short term. Store at -20°C for long term. Avoid repeated freezing-thawing cycles.

BACKGROUND

Introduction	H. pylori is a common type of bacteria that grows in the digestive tract and tends to attack the stomach lining. This bacteria can change the environment around it and reduce the acidity so it can survive more easily. The spiral shape of H. pylori allows it to penetrate the stomach lining, where it's protected by mucus and the body's immune cells can't reach it. This can lead to
--------------	---

stomach problems.

Keywords

H. pylori; Helicobacter pylori
