



Duck Anti-Human ALB Polyclonal Antibody (CABT-CS498)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	ALB
Target	ALB
Isotype	IgY (Δ Fc)
Source/Host	Duck
Species Reactivity	Human
Purification	>95%
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Store at 2-8°C for short term. Store at -20°C for long term. Avoid repeated freezing-thawing cycles.

BACKGROUND

Introduction	<p>Albumin refers generally to any protein that is water soluble, which is moderately soluble in concentrated salt solutions, and experiences heat denaturation. They are commonly found in blood plasma, and are unique to other plasma proteins in that they are not glycosylated. Substances containing albumin, such as egg white, are called albuminoids.</p>
--------------	--

KeywordsALB; albumin
