



Duck Anti-Human AFP Polyclonal Antibody (CABT-CS499)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	AFP
Target	AFP
Isotype	IgY (Δ Fc)
Source/Host	Duck
Species Reactivity	Human
Purification	>95%
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Store at 2-8°C for short term. Store at -20°C for long term. Avoid repeated freezing-thawing cycles.

BACKGROUND

Introduction	Alpha-fetoprotein (AFP) is a protein produced in the liver of a developing fetus. During a baby's development, some AFP passes through the placenta and into the mother's blood. An AFP test measures the level of AFP in pregnant women during the second trimester of pregnancy.
Keywords	AFP; alpha-fetoprotein; α -fetoprotein; FETA; HPAFP