



Duck Anti-Human B2M Polyclonal Antibody (CABT-CS500)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	B2M
Target	B2M
Isotype	IgY (ΔFc)
Source/Host	Duck
Species Reactivity	Human
Purification	>97%
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Store at 2-8°C for short term. Store at -20°C for long term. Avoid repeated freezing-thawing cycles.

BACKGROUND

Introduction	Beta 2-Microglobulin is a low molecular weight protein with sequence homology to immunoglobulins. As a portion of the HLA complex this protein is an important cell-surface structure. Under normal conditions beta 2-microglobulin is synthesized and shed by many cells, particularly lymphocytes, and is detectable in the circulation of normal individuals.
--------------	--

Keywords

B2M; beta-2-microglobulin; B2MG; B2Microglobulin; β 2-Microglobulin; Thymotaxin
