



Duck Anti-Human Ferritin Polyclonal Antibody (CABT-CS503)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Ferritin
Target	Ferritin
Isotype	IgY (Δ Fc)
Source/Host	Duck
Species Reactivity	Human
Purification	>98%
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Store at 2-8°C for short term. Store at -20°C for long term. Avoid repeated freezing-thawing cycles.

BACKGROUND

Introduction	Ferritin is a universal intracellular protein that stores iron and releases it in a controlled fashion. The protein is produced by almost all living organisms, including archaea, bacteria, algae, higher plants, and animals. It is the primary intracellular iron-storage protein in both prokaryotes and eukaryotes, keeping iron in a soluble and non-toxic form. In humans, it acts as a buffer
--------------	---

against iron deficiency and iron overload.

Keywords

Ferritin
