



# Duck Anti-Human IgG Polyclonal Antibody (CABT-CS506)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	IgG
Target	IgG
Isotype	IgY ( $\Delta$ Fc)
Source/Host	Duck
Species Reactivity	Human
Purification	>95%
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Store at 2-8°C for short term. Store at -20°C for long term. Avoid repeated freezing-thawing cycles.

## BACKGROUND

Introduction	Immunoglobulin G (IgG) is a type of antibody. Representing approximately 75% of serum antibodies in humans, IgG is the most common type of antibody found in blood circulation. IgG molecules are created and released by plasma B cells. Each IgG antibody has two paratopes.
Keywords	IgG; Immunoglobulin G;