



Duck Anti-Krait venom Polyclonal Antibody (CABT-CS512)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Krait venom
Target	Krait venom
Isotype	IgY (Δ Fc)
Source/Host	Duck
Species Reactivity	B. multicinctus
Purification	>98%
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Store at 2-8°C for short term. Store at -20°C for long term. Avoid repeated freezing-thawing cycles.

BACKGROUND

Introduction	Krait venom contains β -bungarotoxins, which are presynaptic neurotoxins with phospholipase A2 activity and considered to be the major cause of paralysis. The pre-synaptic action is irreversible and is the reason that once paralysis develops it is not reversed with antivenom
Keywords	Krait venom; Krait