



# Duck Anti-Rabbit IgG Polyclonal Antibody (CABT-CS516)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                                                                                      |
|--------------------|------------------------------------------------------------------------------------------------------|
| Specificity        | IgG                                                                                                  |
| Target             | IgG                                                                                                  |
| Isotype            | IgY ( $\Delta$ Fc)                                                                                   |
| Source/Host        | Duck                                                                                                 |
| Species Reactivity | Rabbit                                                                                               |
| Purification       | >97%                                                                                                 |
| Conjugate          | unconjugated                                                                                         |
| Applications       | ELISA                                                                                                |
| Format             | Liquid                                                                                               |
| Preservative       | None                                                                                                 |
| Storage            | Store at 2-8°C for short term. Store at -20°C for long term. Avoid repeated freezing-thawing cycles. |

## BACKGROUND

|              |                                                                                                                                                                                                                                                                                |
|--------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Introduction | Immunoglobulin G (IgG) is a type of antibody. Representing approximately 75% of serum antibodies in humans, IgG is the most common type of antibody found in blood circulation. IgG molecules are created and released by plasma B cells. Each IgG antibody has two paratopes. |
| Keywords     | IgG; Immunoglobulin G;                                                                                                                                                                                                                                                         |