



Duck Anti-Phenobarbital Polyclonal Antibody (CABT-CS519)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Phenobarbital
Target	Phenobarbital
Isotype	IgY (ΔFc)
Source/Host	Duck
Species Reactivity	N/A
Purification	>97%
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Store at 2-8°C for short term. Store at -20°C for long term. Avoid repeated freezing-thawing cycles.

BACKGROUND

Introduction	Phenobarbital is a prescription medicine used to treat and prevent the symptoms of seizures, sedation, hypnotics, Insomina and Status Epilepticus. Phenobarbital may be used alone or with other medications. Phenobarbital belongs to a class of drugs called Anticonvulsants, Barbiturates.
--------------	---

KeywordsPhenobarbital
