



# Duck Anti-S. Mutans (ATCC25175) Polyclonal Antibody (CABT-CS523)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	S. Mutans
Target	S. Mutans
Isotype	IgY ( $\Delta$ Fc)
Source/Host	Duck
Species Reactivity	S. Mutans
Purification	>95%
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Preservative	None
Storage	Store at 2-8°C for short term. Store at -20°C for long term. Avoid repeated freezing-thawing cycles.

## BACKGROUND

Introduction	Streptococcus mutans is a facultatively anaerobic, gram-positive coccus (round bacterium) commonly found in the human oral cavity and is a significant contributor to tooth decay. It is part of the "streptococci" (plural, non-italic lowercase), an informal general name for all species in the genus Streptococcus.
--------------	--

**Keywords**

S. Mutans; Streptococcus mutans; streptococci

---