



Human Anti-Human CD3 & CD20 (Mosunetuzumab) Monoclonal antibody, clone Mosunetuzumab (CABT-CS550)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	CD3 & CD20
Target	CD3 & CD20
Isotype	IgG1
Source/Host	Human
Species Reactivity	Human
Clone	Mosunetuzumab
Purification	Protein A
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4. Contains no stabilizers or preservatives
Preservative	None
Storage	2 weeks, 2-8°C under sterile conditions after reconstitution. Avoid repeated freeze-thaw. -80°C for a long-term storage.

BACKGROUND

Introduction	The CD3 subunit complex which is crucial in transducing antigen-recognition signals into the cytoplasm of T cells and in regulating the cell surface expression of the TCR complex. T cell activation through the antigen receptor (TCR) involves the cytoplasmic tails of the CD3 subunits CD3 gamma, CD3 delta, CD3 epsilon and CD3 zeta.
Keywords	Cluster of differentiation 3; CD3