



# Human Anti-Human Claudin 18.2 (Zolbetuximab) Monoclonal antibody, clone Zolbetuximab (CABT-CS565)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	CLDN18
<b>Target</b>	CLDN18
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Human
<b>Species Reactivity</b>	Human
<b>Clone</b>	Zolbetuximab
<b>Purification</b>	Protein A
<b>Conjugate</b>	unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	PBS, pH 7.4. Contains no stabilizers or preservatives
<b>Preservative</b>	None
<b>Storage</b>	2 weeks, 2-8°C under sterile conditions after reconstitution. Avoid repeated freeze-thaw. -80°C for a long-term storage.

# BACKGROUND

Introduction	Claudin 18.2 plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity.
Keywords	Claudin-18; CLDN18; CLDN18; claudin 18; SFTA5; SFTPJ