



Human Anti-Human FcRn (Batoclimab)

Monoclonal antibody, clone Batoclimab (CABT-CS569)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	FcRn
Target	FcRn
Isotype	IgG1
Source/Host	Human
Species Reactivity	Human
Clone	Batoclimab
Purification	Protein A
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4. Contains no stabilizers or preservatives
Preservative	None
Storage	2 weeks, 2-8°C under sterile conditions after reconstitution. Avoid repeated freeze-thaw. -80°C for a long-term storage.

BACKGROUND

Introduction	FCGRT binds to the Fc region of monomeric immunoglobulins gamma. It mediates the uptake of IgG from milk. It plays a possible role in transfer of immunoglobulin G from mother to fetus.
Keywords	FcRn; FCGRT; alpha-chain