



# Human Anti-Human Nectin-4 (Enfortumab)

## Monoclonal antibody, clone Enfortumab

### (CABT-CS591)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	PVRL4
<b>Target</b>	PVRL4
<b>Isotype</b>	IgG4
<b>Source/Host</b>	Human
<b>Species Reactivity</b>	Human
<b>Clone</b>	Enfortumab
<b>Purification</b>	Protein A
<b>Conjugate</b>	unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	PBS, pH 7.4. Contains no stabilizers or preservatives
<b>Preservative</b>	None
<b>Storage</b>	2 weeks, 2-8°C under sterile conditions after reconstitution. Avoid repeated freeze-thaw. -80°C for a long-term storage.

# BACKGROUND

Introduction	Poliovirus receptor-like proteins (PVRLs), such as PVRL4, are adhesion receptors of the immunoglobulin superfamily and function in cell-cell adhesion.
Keywords	PVRL4; poliovirus receptor-related 4; Nectin-4; Nectin 4