



Human Anti-Human TIM3 (Sabatolimab) Monoclonal antibody, clone Sabatolimab (CABT-CS601)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	HAVCR2
Target	HAVCR2
Isotype	IgG1
Source/Host	Human
Species Reactivity	Human
Clone	Sabatolimab
Purification	Protein A
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4. Contains no stabilizers or preservatives
Preservative	None
Storage	2 weeks, 2-8°C under sterile conditions after reconstitution. Avoid repeated freeze-thaw. -80°C for a long-term storage.

BACKGROUND

Introduction	HAVCR2 /TIM-3 is member of TIM immunoregulatory proteins family which is encoded by gene on mouse chromosome 11B1.1 and on human chromosome 5q33.2. This chromosomal region has been repeatedly linked with asthma, allergy and autoimmunity
Keywords	HAVCR2; TIM3; TIMD3; FLJ14428; KIM3; TIM-3