



Rabbit Anti-IAV H5N2 (A/American green-winged teal/California/HKWF609/07) HA Polyclonal Antibody (CABT-CS710)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	H5N2 Hemagglutinin/HA
Target	H5N2 HA
Immunogen	Recombinant H5N2 (A/American green-winged teal/California/HKWF609/07) HA / Hemagglutinin protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	IAV
Purification	Protein A & Antigen Affinity
Conjugate	unconjugated
Applications	IF, HI
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS
Preservative	None
Storage	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -

BACKGROUND

Introduction

Influenza hemagglutinin (HA) is a homotrimeric glycoprotein found on the surface of influenza viruses and is integral to its infectivity. HA is a Class I Fusion Protein, having multifunctional activity as both an attachment factor and membrane fusion protein. Therefore, HA is responsible for binding Influenza virus to sialic acid on the surface of target cells, such as cells in the upper respiratory tract or erythrocytes, causing as a result the internalization of the virus. Secondly, HA is responsible for the fusion of the viral envelope with the late endosomal membrane once exposed to low pH (5.0-5.5).

Keywords

H5N2 HA; IAV; IAV H5N2; IAV H5N2 HA; H5N2; Influenza A haemagglutinin H5; H5N2 haemagglutinin
