



Rabbit Anti-IAV H1N1 (A/New Caledonia/20/99) HA1 Polyclonal Antibody (CABT-CS792)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts with human and avian H1(H1N1). Cross-reactivity to other HAs not tested.
Target	H1N1 HA
Immunogen	in vivo expressed hemagglutinin (amino acid 1-343) (H1N1) (A/New Caledonia/20/99) (GenBank No. ABF21272)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	IAV
Purification	Immunoaffinity chromatography
Conjugate	unconjugated
Applications	WB, ELISA
Format	Liquid
Concentration	2 mg/mL
Size	100 µg
Buffer	PBS with less than 0.1% gelatin and 0.1% sodium azide
Preservative	0.1% sodium azide
Storage	Store at -20°C; Do not freeze and thaw. Stable for 6-months from the date of shipment when kept at 4°C. Nonhazardous. No MSDS required.

BACKGROUND

Introduction Influenza hemagglutinin (HA) is a homotrimeric glycoprotein found on the surface of influenza viruses and is integral to its infectivity. HA is a Class I Fusion Protein, having multifunctional activity as both an attachment factor and membrane fusion protein. Therefore, HA is responsible for binding Influenza virus to sialic acid on the surface of target cells, such as cells in the upper respiratory tract or erythrocytes, causing as a result the internalization of the virus. Secondly, HA is responsible for the fusion of the viral envelope with the late endosomal membrane once exposed to low pH (5.0-5.5).

Keywords H1N1 HA; IAV; IAV H1N1; IAV H1N1 HA; H1N1; Influenza A haemagglutinin H1; H1N1 haemagglutinin