



Mouse Anti-MPXV M1R Monoclonal Antibody, Clone 1395 (CABT-CS822)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	MPXV M1R
Target	M1R
Isotype	IgG2a, kappa
Source/Host	Mouse
Species Reactivity	MPXV
Clone	1395
Conjugate	unconjugated
Applications	ELISA, Neut
Format	Liquid
Size	100 µg
Buffer	PBS buffer, pH 7.4
Preservative	None
Storage	4°C for short term (1 week), store at -20°C to -80°C for long term(1 year); Avoid repeated freeze-thaw cycles
Ship	Dry ice

BACKGROUND

Introduction

Monkeypox virus (MPV) is a double-stranded DNA, zoonotic virus and a species of the genus Orthopoxvirus in the family Poxviridae. It is a disease that has similar characteristics to smallpox but with a milder rash and a lower death rate. The virus causes disease in both humans and animals.

Keywords

MPXV; Monkeypox Virus; Monkeypox; MPV
