



Mouse Anti-MNV Polyclonal antibody (Positive Control) (CABT-CS945)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	MNV
Target	MNV
Immunogen	MNV
Isotype	IgG
Source/Host	Mouse
Species Reactivity	MNV
Conjugate	unconjugated
Applications	ELISA
Format	Liquid
Size	1 ml
Storage	Store at 4°C.

BACKGROUND

Introduction	Murine Norovirus (MNV) is a Calicivirus, and unlike the human form which is thought to cause 50% of all viral GI disease worldwide, the mouse form, discovered in 2003, is nonpathogenic (except in very specific strains of immunocompromised mice) and has a prevalence of 58-64% in US populations.
Keywords	Mouse Norovirus; MNV; Norovirus