



# Human Anti-HBeAg Polyclonal antibody (CABT-CS950)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	This polyclonal antibody react with Hepatitis B Virus e antigen
<b>Target</b>	HBeAg
<b>Isotype</b>	IgG
<b>Source/Host</b>	Human
<b>Species Reactivity</b>	HBV
<b>Conjugate</b>	unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M PBS(pH7.4)+ 20%Glycerol
<b>Preservative</b>	0.09% Sodium azide
<b>Storage</b>	Store at 2-8°C up to 21 days, -20°C for long term. Avoid multiple freeze-thaw cycles.

## BACKGROUND

<b>Introduction</b>	HBeAg is a hepatitis B viral protein, produced by the HBcAg reading frame. It is an indicator of active viral replication; this means the person infected with Hepatitis B can likely transmit the virus on to another person (i.e. the person is infectious). HBeAg is considered a marker for cccDNA replication.
---------------------	---

**Keywords**

HBV; Hepatitis B Virus; HBeAg; Hepatitis B Virus e antigen

---