



Magic™ Mouse Anti-Fungal Enolase

Monoclonal antibody, clone K31L (CABT-CS973)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Fungal Enolase
Target	Fungal Enolase
Immunogen	Enolase recombinant protein
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	C. albicans, C. parapsilosis, C. tropicales and C. glabrata, S. cerevisiae.
Clone	K31L
Purification	Protein A
Conjugate	unconjugated
Applications	ELISA (Cap)
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.2
Preservative	0.1% sodium azide
Storage	Store at 4°C

BACKGROUND

Introduction	The glycolytic enzyme enolase is one of the most abundant proteins expressed in fungi and has been shown to be an immunodominant cell-wall-associated antigen of the pathogenic fungus, Candida albicans.
Keywords	Fungal Enolase; Fungal; Enolase; fungi; Candida; Candida albicans; C. albicans