



Magic™ Mouse Anti-Fungal Enolase

Monoclonal antibody, clone I78I (CABT-CS975)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Fungal Enolase
Target	Fungal Enolase
Immunogen	Enolase recombinant protein
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	C. albicans, C. parapsilosis, C. tropicales and C. glabrata, S. cerevisiae.
Clone	I78I
Purification	Protein A
Conjugate	unconjugated
Applications	ELISA (Det)
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.2
Preservative	0.1% sodium azide
Storage	Store at 4°C

BACKGROUND

Introduction The glycolytic enzyme enolase is one of the most abundant proteins expressed in fungi and has been shown to be an immunodominant cell-wall-associated antigen of the pathogenic fungus, *Candida albicans*.

Keywords Fungal Enolase; Fungal; Enolase; fungi; *Candida*; *Candida albicans*; *C. albicans*
