



Rabbit Anti-Human UBE2Q1 polyclonal Antibody (CABT-JX107)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Binds to endogenous Ubiquitin-Conjugating Enzyme E2Q1.
Target	UBE2Q1
Immunogen	Synthesized peptide derived from part region of human protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	The Antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Conjugate	Unconjugated
Applications	WB, ELISA WB 1:500-2000; ELISA 1:5000-20000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	1 mg/mL
Size	20 µg, 50 µg, 100 µg, 200 µg
Buffer	PBS, 50% Glycerol and 0.02% Sodium Azide.
Preservative	0.02% Sodium Azide

Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
----------------	--

Ship	Wet ice
-------------	---------

GENE INFORMATION

Entrez Gene ID	55585
-----------------------	-----------------------

UniProt ID	Q7Z7E8
-------------------	------------------------