



# Rabbit Anti-Beta lactoglobulin Polyclonal Antibody (CABT-JX111)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	LGB
<b>Immunogen</b>	Bovine milk Beta-lactoglobulin (Native Protein)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Bovine
<b>Purification</b>	Protein G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, WB Recommended dilution: WB:1:1000-1:5000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	50 µg, 100 µg
<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Ship

Wet ice

---

## GENE INFORMATION

UniProt ID

[P02754](#)

---