



Quail Anti-Influenza Hemagglutinin-HA peptide polyclonal Antibody (CABT-JX112)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This Antibody recognizes a single band with the predicted molecular weight in Western blot analysis.
Immunogen	Purified recombinant influenza hemagglutinin-HA peptide (YPYDVPDYA).
Isotype	IgY
Source/Host	Quail
Species Reactivity	Influenza virus
Purification	Affinity purified using the immunogen immobilized on a solid support
Conjugate	Unconjugated
Applications	WB, ELISA WB: 1:500 – 1:10,000, ELISA: 1:1,000 – 1:20,000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	1.0 mg/ml
Size	100 µL
Buffer	PBS, pH7.4 with 0.05% Sodium Azide
Preservative	0.05% Sodium Azide

Storage	Store at -20°C. Aliquot to avoid repeated freeze-thaw cycles
Ship	Wet ice