



## Goat Anti-Quail IgY polyclonal Antibody (CABT-JX117)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Target</b>	IgY
<b>Immunogen</b>	Quail IgY (Whole molecule)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	Quail
<b>Purification</b>	Affinity purified using the immunogen immobilized on a solid support
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IA Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Liquid
<b>Concentration</b>	2.0 mg/ml
<b>Size</b>	500 µL
<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C. Aliquot to avoid repeated freeze-thaw cycles

Ship

Wet ice

---