



Rabbit Anti-Rat PTH (aa32-115) Polyclonal Antibody (CABT-JX122)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	The Antibody is a rabbit polyclonal Antibody raised against PTH.
Target	Parathyroid Hormone
Immunogen	Recombinant PTH (Ala32-Gln115) expressed in E. coli.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Antigen-specific affinity chromatography followed by Protein A affinity chromatography
Conjugate	Unconjugated
Applications	WB, ELISA, ICC ELISA 0.05-2ug/ml; ICC 5-20ug/ml; WB 0.5-2ug/ml Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	1.0 mg/ml; 0.5 mg/ml
Size	100 µL, 200 µL, 500 µL, 1 mL
Buffer	0.01 M PBS, pH 7.4, 0.05% ProClin™ 300, 50% Glycerol
Preservative	None

Storage	Avoid repeated freeze/thaw cycles. Store at 4°C for frequent use. Aliquot and store at -20°C for 12 months.
----------------	---

Ship	Wet ice
-------------	---------

GENE INFORMATION

UniProt ID	P01270
-------------------	------------------------
