



Mouse Anti-Canine Haptoglobin Monoclonal Antibody, Clone 10B23 (CABT-JX53009)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Haptoglobin
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Canine
Clone	10B23
Purification	Protein A purified
Conjugate	Unconjugated
Applications	LFIA, ELISA (Cap) We recommend the following as antibody pair (Capture - Detection): CABT-JX53009-CABT-JX53010 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS with 0.1% Sodium Azide
Preservative	0.1% Sodium Azide

Storage	Store at 2-8°C, do not need to freeze.
----------------	--

Ship	Wet ice
-------------	---------
