



# Mouse Anti-Feline Haptoglobin Monoclonal Antibody, Clone 10D22 (CABT-JX53012)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	Haptoglobin
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Feline
<b>Clone</b>	10D22
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	LIA, LFIA, ELISA (Det) We recommend the following as antibody pair (Capture - Detection): CABT-JX53011-CABT-JX53012 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	1 mg
<b>Buffer</b>	PBS with 0.1% Sodium Azide
<b>Preservative</b>	0.1% Sodium Azide

<b>Storage</b>	Store at 2-8°C, do not need to freeze.
----------------	--

---

<b>Ship</b>	Wet ice
-------------	---------

---