



Mouse Anti-Canine Hemoglobin Monoclonal Antibody, Clone 2B4 (CABT-JX53014)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Hemoglobin
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Canine
Clone	2B4
Purification	Protein A purified
Conjugate	Unconjugated
Applications	<p>LFIA (Det)</p> <p>We recommend the following as antibody pair (Capture - Detection): CABT-JX53013-CABT-JX53015</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS with 0.1% Sodium Azide
Preservative	0.1% Sodium Azide

Storage	Store at 2-8°C, do not need to freeze.
----------------	--

Ship	Wet ice
-------------	---------
