



Mouse Anti-Feline pancreatic lipase Monoclonal Antibody, Clone 6H3 (CABT- JX53017)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	PNLIP
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Feline
Clone	6H3
Purification	Protein A purified
Conjugate	Unconjugated
Applications	LFIA, CLIA, ELISA (Cap) We recommend the following as antibody pair (Capture - Detection): CABT-JX53017-CABT-JX53018 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS with 0.1% Sodium Azide

Preservative	0.1% Sodium Azide
Storage	Store at 2-8°C, do not need to freeze.
Ship	Wet ice